

Hemoglobin analysis by tandem mass spectrometry

Christopher A. Haynes, PhD
Service Fellow, CDC


National Conversation: Tandem Mass Spectrometry in Newborn Screening
February 5-6, 2015

OVERVIEW

- Traditional hemoglobin analyses: IEF, HPLC-UV, DNA
- FIA-ESI-MS/MS hemoglobin analysis
- HPLC-ESI-MS/MS hemoglobin analysis
- Future directions

Hemoglobin analysis by isoelectric focusing (IEF)

Identification basis: co-migration



1: β thalassemia trait

2: Normal adult

3: Homozygous β^+ thalassemia

4: Hb A / Hb Korle Bu (D73N)

5: Hb A / Hb D-Punjab (E121Q)

6: Hb A / Hb O-Arab (E121K)

7: Normal newborn

8: Hb A / Hb E (E26K)


9: Hb A / Hb C (E6K)

10: Hb A / Hb S (E6V)

11: Hb S / Hb C

Hemoglobin analysis by HPLC-UV

Identification basis: retention time



Hemoglobin analysis by DNA methodologies

Identification basis: PCR/ microarray

PCR-based methods:

Sanger sequencing of HBA1, HBA2, and HBB

("gold standard"-potentially identifies any mutation)

ASO (Allele-specific oligonucleotide hybridization)

ARMS (Amplification refractory mutation system)

PCR-RFLP (Restriction fragment length polymorphism)

Gap-PCR (PCR amplification across breakpoints)

Microarray-based methods:

Nanochip® (Nanogen, San Diego, CA, USA)

APEX (Arrayed primer extension)

SBE (Single-base extension)

MLPA (Multiple ligation probe amplification)

Hemoglobin analysis by different methodologies

Analysis	Advantages	Disadvantages
IEF	High throughput Inexpensive	Presumptive diagnosis (co-migration)
HPLC-UV	High throughput Inexpensive	Presumptive diagnosis (elution time)
DNA (Microarray, exon seq., Whole-gene seq., etc.)	Definitive diagnosis (amplicon / sequence)	Moderate throughput High equipment cost High reagent cost Expert analyst
ESI-MS/MS (Including trypsin digestion)	Definitive diagnosis (elution time w/ HPLC, precursor <i>m/z</i> , product ion scan) Low reagent cost	Moderate throughput High equipment cost Expert analyst

Hemoglobin: Monomer amino acid sequence

β -globin and α -globin mutants

β MVHLTP**E**EKSAVTALWGKVNVDVGG**E**ALGRLLVYPWTQRFFESFGDLSTPDA
 VMGNPKVKAHGKKVLGAFSDGLAHLNLDKGTFTLSELHCDKLHVDPENFRLL
 GNVLVCVLAHHFGK**E**FTPPVQAAYQKVVAGVANALAHKYH

Globin	Name	HGVS* traditional	HGVS* recommended
Beta	HbS	(c.20A>T,p.E6V)	(c.20A>T,p.E7V)
Beta	HbC	(c.19G>A,p.E6K)	(c.19G>A,p.E7K)
Beta	HbE	(c.79G>A,p.E26K)	(c.79G>A,p.E27K)
Beta	HbD-Los Angeles	(c.364G>C,p.E121Q)	(c.364G>C,p.E122Q)
Beta	HbO-Arab	(c.364G>A,p.E121K)	(c.364G>A,p.E122K)
Alpha	HbG-Philadelphia	(c.[207C>G or 207C>A],p.N68K)	(c.[207C>G or 207C>A],p.N69K)

α MVLSPADKTNVKAAWGKVGAAHAGEYGAELERMFLSFPTTKTYFPHFDLS
 HGSAQVKGHGKKVADALT**N**VAHAVDDMPNALSALSIDLHAHKLRVDPVNF
 KLLSHCLLVTLAAHLPAEFTPAVHASLDKFLASVSTVLTSKYR

Tryptic peptides of β -globin

Trypsin cleaves proteins C-terminal to Arg (R) and Lys (K)



VHLTPEEK

β T1-Met

m/z 952.5*

SAVTALWGK

β T2

m/z 932.5

VNVDEVGGEALGR

β T3

m/z 1314.7

FFESFGDLSTPDAVMGNPK β T5

m/z 2058.9

* m/z values of peptides predicted with Protein Prospector: <http://prospector.ucsf.edu/prospector/mshome.htm>

Tryptic peptides of β -globin mutants

Trypsin cleaves proteins C-terminal to Arg (R) and Lys (K)



VHLTP**VEK**

ST1-Met

m/z 922.5

SAVTALW**GK**

β T2

m/z 932.5

VNVDEV**GK**

ET3

m/z 916.5

FFESFGDLSTPD**AVMGNPK** β T5

m/z 2058.9

* m/z values of peptides predicted with Protein Prospector: <http://prospector.ucsf.edu/prospector/mshome.htm>

Newborn Blood Spot Screening for Sickle Cell Disease by
Using Tandem Mass Spectrometry:
Implementation of a Protocol to Identify Only the Disease
States of Sickle Cell Disease

Stuart J. Moat,^{1*} Derek Rees,¹ Lawrence King,² Adeboye Ifederu,³ Katie Harvey,³ Kate Hall,⁴ Geoff Lloyd,¹
Christine Morell,⁵ and Sharon Hillier⁶

Clin Chem. 2014 Feb ;60(2):373-80

Goal: Identify only sickle-cell disease (HbSS), not carriers (HbAS) or variants

Tryptic digest of DBS punches: SpOtOn kit (Clinical Diagnostics) with labeled peptide

Flow injection analysis: CH₃CN / H₂O (50:50) with 0.025% HCOOH, 2 minutes


MS/MS: Single-reaction monitoring of selected peptides (most are 2⁺ precursors)

Data analysis: ratios of peptide signals

e.g. S:A = (461.7 → 472.2) / (476.7 → 502.2)

Hemoglobin analysis by FIA-ESI-MS/MS (Wales)

Clin Chem. 2014 Feb ;60(2):373-80



Hemoglobin analysis by HPLC-ESI-MS/MS at CDC

- Punch (3.2 mm) DBS; controls are HbA_S adult and HbF_A umbilical cord blood
- Trypsin digestion of punches in 96-well plate (50°C for 30 min)
- Digestion solution contains stable isotope labeled peptide
- Stop digestion with 3 μL of HCOOH per well
- HPLC: Jupiter Protea C12 column (2 x 50 mm), A = 1:99 CH₃CN/H₂O 0.1% HCOOH, B = CH₃CN 0.1% HCOOH, 10 min linear gradient from 0 to 35% B
- MS/MS: Thermo LTQ-XL ion trap, positive ion mode, single-reaction monitoring

α-globin: T1, T10

β-globin: T1, T3, T12

γ-globin: T3, T13

Variants: ST1, CT1, ET3, DT12, OT12, GT10

Internal standard: cut and un-cut (VHLTP[+6]EEKSAVTAL)


Cost: \$0.30 per sample

Product ion scan: m/z 952.5 @6.2 minutes

8.60 E4 cps

β T1-Met

$[M + H - H_2O]^+$
934.5




Blue arrows indicate loss or gain of 18 Da (water)

Hemoglobin analysis by HPLC-ESI-MS/MS at CDC

Variant identification: SRM signal at appropriate elution time
(HbS, HbC, HbE, HbD-Los Angeles, HbO-Arab, HbG-Philly)

Beta-thalassemia identification: quantitation of β T1 peptide

$(952.5 \rightarrow 502.3) / (958.5 \rightarrow 508.3) \times [IS, \mu M] \times \text{Dilution Factor}$



NBS DBS: $944 \pm 399 \mu M$

HbE / β^0 : $511 \mu M$

HbS / β^+ : $85 \mu M$


HbS / β^+ : $134 \mu M$

Analysis of DBS stored under high humidity

B.W. Adam, C.A. Haynes, D.L. Chafin, V.R. De Jesus, Clin Chim Acta, 429 (2014) 59-60.

Challenge: Storage of DBS without desiccant under high humidity conditions

Solution: Analysis of hemoglobin by HPLC-ESI-MS/MS



Hemoglobin analysis by ESI-MS/MS

- Higher specificity than IEF or HPLC-UV analyses
- Faster and less expensive than DNA sequencing
- Simple extraction of DBS punches with trypsin digestion
- Can use FIA-ESI-MS/MS (Wales), ion ratios to identify HbSS
- Can use HPLC-ESI-MS/MS (CDC), quantitation to identify thalassemias
- Suitable for specimens compromised by high humidity

Acknowledgements

Stephanie Guerra
Jessalyn Fontana
David Foreman
Candice Hovell
Victor R. De Jesus
Carla Cuthbert

For more information please contact Centers for Disease Control and Prevention

1600 Clifton Road NE, Atlanta, GA 30333

Telephone: 1-800-CDC-INFO (232-4636)/TTY: 1-888-232-6348

E-mail: cdcinfo@cdc.gov Web: <http://www.cdc.gov>

The findings and conclusions in this report are those of the authors and do not necessarily represent the official position of the Centers for Disease Control and Prevention.